

USER MANUAL

* Images and diagrams referenced in this manual are for illustrative purposes only, The actual appearance may differ.

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

Bullet massage head

Operating Instructions

On/Off: To turn On/Off with a long press on the power button for 1.5s

Intensity level switching: The default level at power-on is Level 1. With some short presses on the power button, the user can switch from Level 1 to Level 6 one by one and then return to Level 1

Timed shutdown: By default, the device will run for 15 minutes before its automatic shutdown, and the running time, 15 minutes, will not be affected by any operations

- Install an appropriate massage head, and switch the intensity level depending on the user's experience. It's recommended to try a low level at the initial use
- Gently press the massage head on the muscle to be relaxed and move it back and forth slowly. If any body part feels particularly sore, apply an appropriate force on it and stay for 15-30 seconds
- Since the body massager generates impact force, the applied pressure should not be too great to cause injury to the muscles and bones
- It's suggested to use the device on a muscular part and stay for max. 60s on the same part
- The recommended duration of use is 15 minutes, which can also be adjusted according to the soreness you feel
- It's recommended to use the device for max. 15 minutes each time, and max. 30 minutes in total

Four replaceable massage heads are available

Ball shaped massage head

Flat massage head

U-shaped massage head

Bullet massage head

- **Ball shaped massage head:** Featuring a large massage area and gentle impact force, and suitable for large muscle groups like the pectoralis major, latissimus dorsi, thigh muscles, and calf muscles
- **Flat massage head:** Featuring hard texture and strong reaction force, and suitable to relax and shape the muscle groups with a certain strength
- **U-shaped massage head:** Used for sore muscles and muscles on both sides of the spine or other special parts
- **Bullet massage head:** Used to reduce pain, suitable for impacting deep tissues of body parts such as palms and soles of the feet

Instructions to Charging/Power Supply

No.	Indicator Status	Status Description
1	Steady green	The battery level is between 25-100% when the device is on
2	Steady orange	The battery level is between 10-25% when the device is on
3	The indicator turns off after 5 orange flashes	The battery level is below 10%
4	The indicator flashes slowly in orange	Charging
5	The indicator flashes slowly in green	Fully charged

Basic Product Specifications

Product name	Body Massager	Rated input power	10W
Model	F1 Pro	Rated input voltage	5V⎓
Battery capacity	1800mAh/13.32Wh		

Routine Maintenance

Keep it away from direct sunlight

Keep it away from high temperature and humidity

Keep it away from smoke and fire

Keep it away from heavy dust

- Keep it away from hard objects and avoid bumping and scratching
- Wipe with a wet tissue or soft cloth for exterior cleaning of the device
- If the body massage will be out of service for a long time, it is recommended to charge it once a month to maintain its sufficient battery level

Precautions for product disposal

Please support the recycling of rechargeable batteries (lithium-ion batteries).

Danger

The battery used in this product is a special rechargeable battery, which should not be used in other products. Do not charge the rechargeable battery that has been taken out.

- Do not overheat the battery by throwing it into the fire
- Do not pierce the battery with nails, impact it by gravity or disassemble it for modification
- Do not connect the positive and negative terminals of the battery directly with metal
- Do not store it together with necklaces, hairpins, etc.
- Do not charge, use and place the battery in a hot place, such as by the fireplace or under the hot sun

All of the above can cause heat, fire, rupture, etc.

Warning

Please keep the removed rechargeable battery out of the reach of infants and children to prevent accidental swallowing.

- In case of accidental swallowing, please seek immediate medical treatment
- When the battery fluid leaks, do not touch the liquid with your bare hands, and it can be processed in the following way.
- When the battery fluid enters the eyes, there is a risk of blindness. Do not rub your eyes directly with your hands, but wash them with clean water immediately and seek medical treatment.
- When your body or clothing is contaminated by the battery electrolyte, it may cause skin inflammation or injury. Please rinse with water. If you still feel discomfort, seek medical treatment immediately.
- Only the product can be disassembled at the time of disposal. Failure to comply may result in fire, electric shock, injury and other accidents. If repairs are required (except for battery replacement), please consult the store where you purchased the product.

The lithium-ion battery used in this product is a valuable resource that can be recycled. Therefore, please take out the lithium-ion battery when the device is disposed of, and hand it over to a professional organization that can recycle rechargeable batteries for disposal.

- Batteries must be removed from the device before it is discarded;
- The device must be powered off when the battery is removed;
- Batteries should be disposed of safely.

Battery removal steps:

1. Use the tool to open the shell;
2. Use scissors to cut off the battery leads separately;
3. Remove the battery;

- * Keep it away from hard objects, avoid bumping and scratching.
- * Wipe the exterior of the device with a soft cloth dipped in an appropriate amount of alcohol to clean it.
- * If the device is not used for a long time, it is recommended to charge it once a month to keep the battery level above 20%.

Safety instructions

EN

- It is not recommended to massage the skin surface directly with this product, so as not to cause damage to the skin. Wear dry and clean clothes or cover the massage site with a clean and dry towel, then gently press and move the massage gun
- Use high speed and high pressure on the soft tissue parts of the body only when there is no pain and discomfort, and do not use high speed and high pressure on the head and any bone parts of the body
- Frequent high-speed or high-pressure application on the same part may cause abrasions, and stop using the massage gun immediately if there is pain and discomfort
- Keep your fingers, hair or any part of the body away from the retractable part of the device to prevent pinch pain and injury
- Do not discard and misuse this product
- Do not make any modifications to this product without authorization
- Do not leave the product unattended when operating or charging it
- This product is intended for adult use only
- People with contraindications should use this product under the permission of a doctor

Note

- Before starting this product, please make sure that there is no obstacle in front of the massage head. If it is blocked by an external force, pause shortly before restarting
- This product is not waterproof. Please do not expose it to water or other liquids
- Do not disassemble the product. If it is disassembled by mistake and works abnormally, please send it back for commissioning. All costs arising therefrom shall be borne by the user
- Continuous work for a long time may cause the surface temperature of the motor to rise. Please operate with caution
- Do not drop or impact the product; otherwise, it may damage, causing abnormal operation. If the product malfunctions due to dropping or impact, please contact the after-sales service in time

Warnings

EN

! Contraindications

- Not for use by cancer patients
- Not for use by patients receiving internal fixation with metal spicula (it is forbidden to use the device directly on the fracture site)
- Not for use by patients with osteoporosis
- Not for use by patients with venous thrombosis and pulsation
- Not for use by patients with heart disease and high blood pressure
- Use with caution by pregnant women (do not use it near the abdomen)
- Use with caution by patients with kidney stones (do not use it near the kidney area)
- Not for use by patients with vasodilation
- Not for use by patients with acute inflammation
- Not for massaging the patella, tibia, metatarsal and other bone parts
- Not for massaging the spine directly
- Not for massaging the cervical spine directly
- Not for massaging the facial muscles directly
- Not for use by patients with medical electronic devices implanted in the body, such as cardiac pacemakers
- Not for use in the states of fasting, a full stomach and drunk

- Do not massage the same part for more than 4-5 minutes
- Please consult a doctor before using for patients with other diseases such as cervical and lumbar spondylosis
- During use or when starting up, do not insert your fingers into the massage head assembly hole to avoid injury from being hit.

♥ Applications

- Relaxation before and after physical training
- Rehabilitation training
- Muscle pain and fascia fatigue relief
- Help with physical training
- Muscular fascia relaxation and soft tissue extensibility improvement

⚠ Danger

- Do not use while charging

⚠ Suggestions:

- Ensure that the main unit buttons are functioning properly before use.
- To prevent pinching injuries, be cautious to avoid placing fingers or other body parts in narrow parts of the device.
- When disposing of this product, please do so in accordance with local environmental protection regulations to prevent environmental pollution.
- Recommended for use by individuals aged 16 and above.

⚠ Please do not use or store this product in the following locations:

- Places exposed to direct sunlight.
- Areas with high temperatures and humidity.
- Locations with the presence of water.
- Dusty environments.
- Near open flames.
- In strong electromagnetic fields.
- Places where the product may be subject to vibrations or shocks.

⚠ Caution:

- If you experience skin irritation when using the body massager, please refrain from using it, and wait two to three days to see if your symptoms ease. If symptoms persist or worsen, please consult a physician.
- If you have eczema, allergies or asthma, you may be more likely to experience a skin irritation or allergy from using the massager.
- * The appliance it must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
- * This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- * Children shall not play with the appliance.
- * Cleaning and user maintenance shall not be made by children without supervision.

CE This device in compliance with the essential requirements and other relevant provisions of Directive (EU)2023/988.
CAUTION : RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE. DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.

EU Compliance Statement:
 SKG Health technology Co.,Ltd. Here by declares that this device is in compliance with the essential requirements and other relevant provisions of the Directive (EU)2023/988. A copy of the EU Declaration of Conformity is available online at <https://www.skg.com/pag-es/skg-declaration-of-conformity>.

Troubleshooting

If the product does not work properly, check for the following problems. If the problem persists, return it to the dealer for service

Problem	Possible Cause	Action
Failure to start	The battery level is low	Charge the product
	The product is not turned on	Press and hold the On/Off button to power on
	The product is still connected to the household power socket	Remove the charging cable from the massager
Failure to charge or the charging light is not On	The charging cable is not connected to the power outlet correctly	Replug the charging cable
	An incorrect charging cable is used	Use a USB-Type-C cable

Old electrical appliances must not be disposed of together with the residual waste, but have to be disposed of separately. The disposal at the communal collecting point via private persons is for free. The owner of old appliances is responsible to bring the appliances to these collecting points or to similar collection points. With this little personal effort, you contribute to recycle valuable raw materials and the treatment of toxic substances.

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SKG Health Technologies Co., Ltd.