
SKG Neck Massager PS700-2

INSTRUCTION MANUAL

SKG

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WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY WHEN CHILDREN ARE PRESENT, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:

READ ALL INSTRUCTIONS BEFORE USING

DANGER – TO REDUCE THE RISK OF ELECTRIC SHOCK:

Always unplug the device from the electrical outlet immediately after using and before cleaning.

- DO NOT reach for a device that has fallen into water. Unplug it immediately.
- DO NOT use the device while bathing or in a shower.
- DO NOT place or store the device where it can fall or be pulled into sink or tub.
- DO NOT place or drop the device in water or other liquid.
- DO NOT plug or unplug the charging cable with wet hands.

WARNING – TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK, FIRE OR INJURY TO PERSONS:

- The device should NEVER be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts or attachments.
- Close supervision is necessary when this device is used by, on, or near children, invalids or disabled persons.
- Children should not use this device.
- Use this device only for its intended use as described in this manual. DO NOT use attachments not recommended by SKG; specifically any attachments not provided with the device.
- NEVER operate this device if it has a damaged power button, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- NEVER drop or insert any object into any opening of the device.
- DO NOT operate the device where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT operate the device under a blanket or pillow. Never cover the device when it is in operation, as excessive heating can occur, causing fire, electric shock or injury to persons.
- To disconnect, turn all controls to the "off" position.
- NEVER operate the device with the air openings blocked. Keep the air openings free of lint, hair, and the like.
- NEVER operate the device on a soft surface such as a bed or couch where the air openings may be blocked.
- Use heated surfaces carefully. DO NOT use over insensitive skin areas or in the presence of poor circulation. The unattended use of the heated device by children or incapacitated persons may be dangerous.
- DO NOT use this massager in close proximity to loose clothing, jewelry, or long hair.

SAVE THESE INSTRUCTIONS

CAUTION—PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE USE.

- This product is not intended for medical use. It is intended only to provide a luxurious massage.
- Consult your doctor prior to using this product if you are pregnant, you have a pacemaker or have any concerns regarding your health.
- This product is not recommended for use by diabetics.

- NEVER leave the device unattended, especially when children are present.
- DO NOT use this product for more than 15 minutes at a time.
- NEVER use this product directly on swollen or inflamed areas, or on skin eruptions.
- DO NOT use this product as a substitute for medical attention.
- NEVER use this product while lying flat or in bed.
- NEVER use this product with a necklace on.
- This product should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- This product should not be used by children or invalids without adult supervision.
- This product is not intended for use in automobiles.
- This product is intended for household use only.

BASIC PRODUCT PARAMETERS

EN

Product Name: Neck Massager

Model: PS700-2

Net Weight: 594g

Packing Size: 268x202x133mm

Battery Capacity: 1400mAh/10.36Wh

Rated Input Voltage: 5V

Rated Input Power: 14W

Operating environment:

Temperature range: 0°C-40°C

Relative Humidity: <80% RH, non-condensing


Atmospheric Pressure: 50Kpa-106Kpa

Storage and transportation environment:

Temperature range: -20°C to +55°C

Relative Humidity range: <95% RH, non-condensing

Atmospheric Pressure: 50Kpa-106Kpa

Note: Images are for reference only. Actual product may vary.

What is in the box?

The PS700-2 device x1
User Manual x1

USB-C charging cable x1
Quick Start Guide x1

Buttons & Controls

• Startup Mode

Upon startup, the device resumes your last used massage mode and heat setting.

• Switching Heat Levels

To change the heat setting, briefly press the "Hot Compress Button". It cycles in the following order:

Medium → High → Off → Low (and repeats).

• Switching Massage Modes

To change the massage mode, briefly press the "Power/Mode Button". It cycles in the following order:

De-stress Mode → Meditation Mode → Relax Mode → Shiatsu Mode (and repeats).

• Auto Shut-off

Massage mode on: Shuts off (indicator lights off) after 15 minutes.

Indicator Lights

Light Color	Light State	Description
Green	Steady On	Battery level $\geq$ 25%
Orange	Steady On	Battery level $<$ 25%
Orange	Flashes, then turns off	Battery level $<$ 9%
Orange	Breathing (slow pulsing)	Charging in progress
Green	Breathing (slow pulsing)	Fully charged

Note: This product does not support charging while in use. If you need normal operation, unplug the charging cable.


Charging

• Before charging, ensure the device is turned off. Use the provided USB-C cable to connect the device to a power adapter for charging.

• A full charge can support approximately 6 complete massage sessions (actual duration may vary depending on the heat level and massage mode selected).

Note: If using a third-party charger, ensure it is certified and in good condition. Upon startup, the device resumes your last used massage mode and heat setting.


Usage Instructions

Step 1:

Wear Correctly


Ensure the massager is properly positioned around your neck. Tip: Remove necklaces and earrings before use. Keep hair away from the massage nodes to prevent entanglement.

Step 2:

Adjust for Fit


Gently adjust the device so the massage nodes fit snugly against your neck.

Step 3:

Start Your Massage

Press and hold the power button, then select your preferred intensity to begin.

Usage Tip

The hand straps are detachable. Pulling them adjusts the fit, allowing you to fine-tune the massage pressure.

The device features dual control via both the mobile app and physical buttons. Through the app, users can access additional massage modes, customize session durations, and intelligently adjust pulse intensity, vibration strength, and heat levels for a personalized massage experience.

SKG Health APP Setup Guide

Step 1: Download the App

Download the "SKG Health" app from the Google Play store or Apple App Store, or scan the provided QR code.

Google Play and the Google Play logo are trademarks of Google Inc.
Apple and the Apple logo are trademarks of Apple Inc.

Step 2: Register/Log In

Open the app and register for a new account or log in to your existing account.

Step 3: Add Your Device

Power on your device and enable Bluetooth on your phone. Then, in the "SKG Health" app, tap the "Add device" button to connect.

Device Maintenance

The following maintenance instructions are important to ensure that your device continues to work as designed. Failure to follow these instructions may cause your device to stop working.

Cleaning the Device

Read the following instructions carefully before cleaning the device to ensure longevity of the device and its components.

1. Visually inspect the device for any obvious signs of debris build-up.
2. Wipe the device with a non-abrasive, lint-free damp cloth or alcohol-free cleansing wipe.
3. Dry the device with a non-abrasive, lint-free cloth.
4. After cleaning, allow the device to dry thoroughly before storing or beginning another treatment. The device is NOT waterproof or machine washable. Do not allow excess liquids or water to enter the product.
5. A properly cleaned device should have no visible signs of debris or moisture.

****Note:** This product is not waterproof. Only clean with a damp cloth or alcohol-free cleansing wipe. Do not submerge the device in water or clean it under running water. Do not allow the device to come into contact with corrosive solutions, which can damage the device's appearance and function.

TROUBLESHOOTING

Refer to the table for solutions to any problems. If the problem persists, please return the product to the dealer for repair.

Condition	Possible Cause	Solution
Won't Power On	Low battery	Charge the device
	The device is currently charging	Unplug the charging cable before use
Won't Charge / No Indicator Light	The charging cable is not securely connected.	1. Unplug the cable. 2. Firmly reconnect it to both the device and the power adapter.

WARNING

No modification of this device is allowed.
Do not modify this device without authorization of the manufacturer.
Any modification requires proper inspection and testing to ensure the device remains safe for use.

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DANGER

Please support the recycling of rechargeable battery (lithium-ion battery).

The battery included in this product is a dedicated rechargeable battery, not intended in other products. Please do not charge the removed rechargeable battery.

- DO NOT throw it into a fire, or heat it up to high temperature.
- DO NOT pierce it with nails, impact it by gravity or disassemble it for modification.
- DO NOT connect the positive and negative poles of the battery directly with metal.
- DO NOT store it together with necklaces, hairpins, etc.
- DO NOT charge, use or place the battery in areas with high temperatures, such as near a fireplace or under direct sunlight.
- The above situations may cause heating, fire and breakage.

WARNING

Keep the removed rechargeable batteries out of reach of children to prevent accidental ingestion.

If a battery is swallowed, consult a doctor immediately.

In the event of a battery fluid leak, avoid touching the fluid and handle according to safety procedures.

- If battery electrolyte enters the eyes, rinse with clean water and seek medical advice without rubbing your eyes.
- If battery electrolyte contacts the skin, rinse the area with water and get medical help if needed.

Only disassemble scrapped products to prevent accidents such as fires or electric shocks. For battery repair, contact the dealer where you purchased the product. This does not include battery replacement.

The lithium-ion battery in this product is recyclable and should not be replaced; surrender it to a professional recycling agency when the product is discarded.

Always turn off the device before removing the battery for disposal.

To safely dispose of the battery, follow these simplified steps:

Use specialized tools to detach the device's casing.

Remove the internal circuit board and motor.

Extract the battery, sever its wires with pliers, and remove it completely.

RF warning statement:

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED STATEMENT

English:

This device complies with Innovation, Science and Economic Development Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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This device was tested for typical body-worn exposure requirements, a minimum separation distance of 0cm must be maintained between the user's body and the handset, including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn devices that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna. This device is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU. All essential radio test suites have been carried out.

1. CAUTION : RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE. DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS

2. The device complies with RF specifications when the device is used in accordance with the body EU Compliance Statement:

SKG Health technology Co.,Ltd. Here by declares that this device is in compliance with essential requirements and other relevant provisions of the Directive 2014/53/EU. A copy of the EU Declaration of Conformity is available online at <https://www.skg.com/pages/skg-declaration-of-conformity>

Old electrical appliances must not be disposed of together with the residual waste, but have to be disposed of separately. The disposal at the communal collecting point via private persons is for free. The owner of old appliances is responsible to bring the appliances to these collecting points or to similar collection points. With this little personal effort, you contribute to recycle valuable raw materials and the treatment of toxic substances.

REGISTER & SUPPORT

EN

Register your SKG product at support@skg.com to maximize support and visit <https://www.skg.com/> for expert tips, tutorials, and FAQs to optimize your device use.

SUCCESS COURIER SL
CALLE RIO TORMES NUM.1, PLANTA 1, DERECHA
OFICINA 3, Fuenlabrada, Madrid, 28947 Spain
succeservice2@hotmail.com

C05_EU_202310214000009

EVATOST CONSULTING LTD
Office 101 32 Threadneedle Street, London,
United Kingdom, EC2R 8AY
contact@evatost.com

UK0017177407401980mk

SKG Health Technologies Co., Ltd.

Address: 23A Floor, Building 3, Zhongke R&D Park, No. 009,
Gaoxin South 1st Road, High-tech Zone Community, Yuehai Street,
Nanshan District, Shenzhen, Guangdong, China.

Website : <https://www.skg.com/>

Made in China